

# Cholinesterase Detection Kit

For Research Use Only.




Figure 1

For technical questions and orders, please contact us at:

1-800-829-3194  
952-888-8788  
952-888-8988 fax  
www.immunochemistry.com

### Cholinesterase Detection Kit Ordering Information:


Kit	Reagent	25-test kit	100-test kit
Cholinesterase	Ph-Fluorescein	catalog# 973	catalog# 974

**FLICA™ Ordering Information:** Use the cholinesterase kit with ICT's red fluorescent FLICA™ Apoptosis Detection Kits for dual-labeling studies. The Poly-Caspases kit (#916 and #917) was used in Figure 2 on page 14:

Caspase	FLICA™ Reagent	25-test kit	100-test kit
Poly-Caspases	SR-VAD-FMK	catalog# 916	catalog# 917
Caspases 3&7	SR-DEVD-FMK	catalog# 931	catalog# 932
Caspase 9	SR-LEHD-FMK	catalog# 960	catalog# 961

**Cover photo, Figure 1:** Localization of cholinesterase with Ph-FI in the nerve-muscle junctions (end-plates) in C57 mouse diaphragm muscle tissue. The diaphragm muscle was dissected and fragments of it incubated in PBS containing 20µM Ph-FI for 1 hour. The muscle was then rinsed for 20 minutes in PBS, stretched on a microscope slide, mounted under a coverslip, and examined by fluorescence microscopy using blue light excitation. Notice the retention of the fluorochrome in the structures with typical features of end-plates (arrow).

No part of this manual may be reprinted without the express written permission of Immunochemistry Technologies, LLC. Publication Version: Cholinesterase Flat Manual #819-030105



www.immunochemistry.com  
Immunochemistry Technologies, LLC  
9401 James Avenue South, Suite 155  
Bloomington, MN 55431  
800-829-3194 / 952-888-8788  
fax 952-888-8988

## Table of Contents

1. Introduction .....	4
Overview of the Cholinesterase Protocol	
2.....	6
3. Contents of the 25-test Cholinesterase Detection Kit #973 .....	6
4. Contents of the 100-test Cholinesterase Detection Kit #974 .....	7
5. Recommended Materials and Equipment .....	7
Instrumentation	
6.....	7
7. Storage and Shelf-Life .....	7
Safety Information	
8.....	8
9. Preparation of 1X Cellular Wash Buffer .....	8
10. Fixative .....	9
11. Reconstitution of the Ph-FI Vial .....	9
12. Storage of the 255X Ph-FI Stock Concentrate for Future Use .....	9
13. Preparation of the 51X Ph-FI Working Solution for Immediate Use ...	10
14. General Procedure .....	10
15. Washing Non-Adherent Cell Populations .....	12
16. Washing Adherent Cell Populations .....	13
17. Sample Microscope Data Figure 2.....	14
18. References .....	15

**For Research Use Only. Not for use in diagnostic procedures.**

## 1. Introduction

Physostigmine and its related compounds have long been known to bind to cholinesterase-type enzymes (acetylcholinesterase and butyrylcholinesterase) leading to their widely accepted use in the study and treatment of cholinergic-neurotransmission-related diseases such as Alzheimer's Disease (AD) (3). As AD progresses, acetylcholinesterase (AChE) activity decreases in some brain regions, while butyrylcholinesterase (BChE) activity increases in other brain regions (4, 5). This may be the result of a relative increase in the numbers of BChE-positive neurons. This increase in BChE activity is greater in the hippocampus and temporal lobe in patients with AD. Also, AChE activity was previously detected in several different apoptotic cell lines and was proposed to be associated with the apoptosome formation during the initial stage of apoptosis (1, 2).

Deposits of amyloid material can be observed many years prior to the actual occurrence of neuronal degeneration and dementia. Both AChE and BChE are associated with amyloid and neurofibrillary tangles in the brains of elderly persons with, and without, relevant cognitive impairment (6). It is therefore important to identify what factors may contribute to the transformation of "benign amyloid" to one capable of causing the classical neurological symptoms of the AD state. BChE is a possible factor in this transformation process. Advanced amyloid plaques in Alzheimer brains have up to 87% BChE reactivity, compared with 20% reactivity in early, benign deposits (7). The presence of BChE activity distinguishes the neuro-toxic plaques from those seen in normal aging. This could indicate that BChE activity plays a role in the transformation of benign amyloid plaques into the forms associated with neuronal degeneration and clinical dementia.

Originally isolated from the Calabar bean vine, *Physostigma venenosum*, physostigmine (a carbamate ester alkaloid) mimics the natural substrate, acetylcholine, thus allowing it to act as a cholinesterase inhibitor (8). Cholinesterase inhibition is facilitated by the ability of the physostigmine inhibitor to target the cholinesterase active site and subsequently carbamylate the reactive serine residue within the active site region of the enzyme (3). The carbamylated enzyme intermediate is much more stable than the acylated-enzyme intermediate which results from its reaction with the natural substrate, acetylcholine. This feature allows us to utilize physostigmine as a cell-permeant cholinesterase-targeting probe when labeled with a green fluorescent tag (fluorescein).

In this assay, fluorescein has been conjugated to a physostigmine analog, eseroline, through a 5-carbon spacer linked to the carbamoyl carbonyl group of the physostigmine side chain. The resulting conjugate (Physostigmine – Fluorescein, Ph-FI) can then be used to detect the activity of cholinesterase enzymes, as well as apoptosis activity in several different cell lines (9).

When the Ph-FI probe reacts with cholinesterase, the carbamoyl-5 carbon-fluorescein remains bound to the serine hydroxyl active site of cholinesterase in the form of a stable intermediate, while simultaneously liberating an eseroline molecule as the hydrolysis product. Eventually, the carbamoyl-5 carbon-fluorescein tag is released from the enzyme active site serine, regenerating the active cholinesterase enzyme form. The bound-cholinesterase form of this complex, following a brief wash step, can be quantitated using 3 different fluorescence detection methods: 96-well microtiter plate fluorometry for quantitation; fluorescence microscopy for qualitative analysis; and flow cytometry for quantitation.

Using a fluorescence plate reader (with **black** microtiter plates), cholinesterase activity can be quantitated as the amount of green fluorescence emitted from bound Ph-FI probes. Cell populations with increasing cholinesterase activity will have a higher RFU intensity than cell populations with less activity.

Viewing cells through a fluorescence microscope, cholinesterase-positive cells will fluoresce green, while negative cells will appear mostly unstained.

Using a flow cytometer, analysis is done using a 15 mW argon ion laser at 488 nm. Fluorescein is measured on the FL1 channel, and a log FL1 (X-axis) versus number of cells (Y-axis) histogram may be generated. On this histogram, there will appear two cell populations represented by two peaks. The majority of the negative cells will occur within the first log decade of the FL1 (X) axis (first peak), whereas the cholinesterase-positive cell population will appear as a separate peak or as a shoulder of the first peak showing increased fluorescence intensity.

The Ph-FI reagent has an optimal excitation range from 488 - 492 nm, and emission range from 515 - 535 nm (the excitation / emission pairs which best approximate this optimal range should be used). Cells labeled with the Ph-FI reagent may be read immediately or preserved for 24 hours using the fixative.

Following the suggested protocols listed here, each 500  $\mu$ L sample of your cell culture (grown up to  $1 \times 10^6$  cells/mL) requires 10  $\mu$ L of the 51X Ph-FI working solution (equal to 2  $\mu$ L of the 255X Ph-FI stock concentrate). Your cells may require more or less of this reagent, and the length of the incubation time may vary based on your protocol.

**For Research Use Only. Not for use in diagnostic procedures.**

## 2. Overview of the Cholinesterase Protocol

Labeling cells with ICT's Cholinesterase Detection Kit can be completed within a few hours. However, the Cholinesterase Detection Kit is used with living cells, which require periodic maintenance and cultivation several days in advance. In addition, once the proper number of cells has been cultivated, time must be allotted for your experimental manipulation of the cells. Therefore, as the 51X Ph-FI working solution must be used immediately, the Ph-FI reagents should be prepared at the end of your experimental process. The following is a quick overview of the Cholinesterase Detection Kit protocol:

1. Culture cells to a density optimal for your specific experiment, but not to exceed  $10^6$  cells/mL.
2. Expose your cells to your experimental conditions.
3. At the same time, culture a non-treated negative control cell population at the same density as the experimental population for every labeling condition. For example, if labeling with Ph-FI and SR-VAD-FMK (from FLICA™ kit #916), make at least 2 populations:
  - a. Treated cells labeled with Ph-FI and SR-VAD-FMK reagents.
  - b. Untreated cells labeled with Ph-FI and SR-VAD-FMK reagents.
4. Dilute the cellular wash buffer in DI H<sub>2</sub>O to 1X (Section 9).
5. Reconstitute the Ph-FI reagent with 50  $\mu$ L DMSO to create the 255X Ph-FI stock concentrate (Section 11).
6. Add 200  $\mu$ L PBS to the reconstituted reagent to create the 51X Ph-FI working solution (Section 13).
7. Add 10  $\mu$ L of the 51X Ph-FI working solution to 500  $\mu$ L of your cells to label them, yielding a final concentration of 20  $\mu$ M Ph-FI at 1X in the cells. If also labeling with SR-FLICA™, add that reagent at the same time.
8. Incubate cells for about 1 hour.
9. Remove the Ph-FI-containing media and wash cells (Section 14, 15, or 16).
10. Add fresh 1X cellular wash buffer, media, or PBS.
11. If desired, fix cells (Section 10).
12. Analyze data via microtiter plate fluorometry, fluorescence microscopy, or flow cytometry.

## 3. Contents of the 25-test Cholinesterase Detection Kit, catalog# 973:

- 1 amber vial of lyophilized Cholinesterase reagent, Ph-Fluorescein (Ph-FI), 25 tests per vial, catalog# 6156.
- 1 bottle of 15 mL 10X Cellular Wash Buffer, catalog# 6164
- 1 bottle of 6 mL Fixative, catalog# 636
- Assay Manual, catalog# 820
- MSDS sheets

**4. Contents of the 100-test Cholinesterase Detection Kit, catalog# 974:**

- 4 amber vials of lyophilized Cholinesterase reagent, Ph-Fluorescein (Ph-FI), 25 tests per vial, catalog# 6156.
- 1 bottle of 60 mL 10X Cellular Wash Buffer, catalog# 6165
- 1 bottle of 6 mL Fixative, catalog# 636
- Assay Manual, catalog# 820
- MSDS sheets

**5. Recommended Materials and Equipment (not all are required):**

- Cultured cells with media
- Reagents to do your experiment
- 15 mL polystyrene centrifuge tube (1 per sample)
- Amber vials or polypropylene tubes for storage of the 255X stock concentrate at  $-20^{\circ}\text{C}$ , if aliquoted
- 150 mL or 600 mL graduated cylinder
- Slides
- Hemocytometer
- Centrifuge at  $<300 \times g$
- $37^{\circ}\text{C}$   $\text{CO}_2$  incubator
- Vortexer
- Pipette(s) capable of dispensing at  $10\mu\text{L}$ ,  $50\mu\text{L}$ ,  $200\mu\text{L}$ ,  $300\mu\text{L}$ ,  $1\text{mL}$
- $\text{DIH}_2\text{O}$ , 135 mL or 540 mL needed
- Phosphate Buffered Saline (PBS) pH 7.4, up to 100 mL needed
- Dimethyl Sulfoxide (DMSO),  $50\mu\text{L}$  or  $200\mu\text{L}$  needed
- Ice or  $4^{\circ}\text{C}$  refrigerator to store cells

**6. Instrumentation (not all are required):**

- 96-well fluorescence plate reader with excitation at 488 nm, emission 520 nm filter pairings, and **black** round or flat bottom 96-well microtiter plates.
- Fluorescence microscope with appropriate filters (excitation 490 nm, emission  $>520$  nm) and slides.
- Flow cytometer equipped with a 15 mW, 488 nm argon excitation laser, with appropriate filters (excitation 490 nm, emission  $>520$  nm).

**7. Storage and Shelf-Life**

- Store the unopened kit (and each unopened component) at  $2^{\circ}\text{C}$  to  $8^{\circ}\text{C}$  until the expiration date.
- Protect the cholinesterase reagent Ph-FI from light at all times.
- Once reconstituted, the 255X stock concentrate should be stored at  $-20^{\circ}\text{C}$  protected from light. This reagent is stable for 6 months and may be thawed twice during that time.
- Once diluted, store the 1X cellular wash buffer at  $2 - 8^{\circ}\text{C}$  for 30 days.

- Additional cellular wash buffer and fixative may be ordered by calling ICT at 1-800-829-3194 or 952-888-8788.

**8. Safety Information**

- Use gloves while handling the Ph-FI reagent, cellular wash buffer, and fixative.
- Dispose of all liquid components down the sink and flush with copious amounts of water. Solid components may be tossed in standard trash bins.
- MSDS sheets are available at [www.immunochemistry.com/MSDS.htm](http://www.immunochemistry.com/MSDS.htm) or by calling 1-800-829-3194 or 952-888-8788.

**9. Preparation of 1X Cellular Wash Buffer**

The cellular wash buffer is supplied as a 10X concentrate which must be diluted to 1X with DI  $\text{H}_2\text{O}$  prior to use.

1. If necessary, gently warm the 10X concentrate to completely dissolve any salt crystals that may have come out of solution. Do not let it boil.
  - a. For the 25-test kit, add the entire bottle (15 mL, catalog# 6164) of 10X cellular wash buffer to 135 mL of DI  $\text{H}_2\text{O}$  (to make 150 mL).
  - b. Or, for the 100-test kit, add the entire bottle (60 mL, catalog# 6165) of 10X cellular wash buffer to 540 mL of DI  $\text{H}_2\text{O}$  (to make 600 mL).
  - c. Or, if not using the entire bottle, dilute the 10X cellular wash buffer 1:10 in DI  $\text{H}_2\text{O}$ . For example, add 10 mL 10X wash buffer to 90 mL DI  $\text{H}_2\text{O}$  (to make 100 mL).
2. Let the solution stir for 5 minutes or until all crystals have dissolved.
3. If not using the 1X cellular wash buffer the same day it was prepared, store it covered at  $2^{\circ} - 8^{\circ}\text{C}$  for 30 days.

The cellular wash buffer contains a small amount of sodium azide which should not affect the cells during the short amount of time they are in the wash buffer. However, if you do not want to use ICT's cellular wash buffer, just use your own fresh cell culture media instead to wash the cells. We do not recommend using PBS to wash the cells. If more cellular wash buffer is needed, please contact ICT at 1-800-829-3194 or 952-888-8788 to order another bottle.

**Warning:** The wash buffer contains sodium azide, which is harmful if swallowed or absorbed through the skin. Sodium azide can react with lead and copper sink drains forming explosive compounds. When disposing of excess wash buffer, flush sink with copious amounts of water; see MSDS for further information.

## 10. Fixative

If the stained cell populations cannot be evaluated immediately upon completion of the staining protocol, cells may be fixed and analyzed up to 24 hours later on a microscope or flow cytometer. The fixative is a formaldehyde solution designed to cross-link cell components and will not interfere with the carboxyfluorescein labeling once the labeling reaction has taken place.

After labeling, add the fixative into the cell solution at a 1:10 ratio. For example, add 100  $\mu\text{L}$  fixative to 900  $\mu\text{L}$  cells. Fixed cells may be stored on ice or at 2-8°C up to 24 hours.

- **Do not use ethanol-based or methanol-based fixatives to preserve the cells - they will inactivate the Ph-FI label.**
- **Never add the fixative until the staining and final wash steps have been completed.**

## 11. Reconstitution of the Ph-FI Vial

The Ph-FI reagent is supplied as a highly concentrated lyophilized powder (the amber vial may appear empty as the reagent is lyophilized onto the insides of the vial). It must first be reconstituted in DMSO, forming a 255X Ph-FI stock concentrate, and then diluted 1:5 in PBS to form a final 51X Ph-FI working solution. For best results, the 51X Ph-FI working solution should be prepared immediately prior to use; however, the reconstituted 255X Ph-FI stock concentrate can be stored at -20°C for future use.

- **The newly reconstituted 255X Ph-FI stock concentrate must be used or frozen immediately after it is prepared and protected from light during handling.**
1. Reconstitute each vial of lyophilized Ph-FI (catalog #6156) with 50  $\mu\text{L}$  DMSO. This yields a 255X stock concentrate. (The 25-test kit contains 1 vial; the 100-test kit contains 4 vials.)
  2. Mix by swirling or tilting the vial, allowing the DMSO to travel around the base of the amber vial until completely dissolved. At room temperature (RT), this reagent should be dissolved within a few minutes.
  3. If immediately using this solution, dilute it to 51X (Section 13).
  4. Or, if using later, aliquot and store it at -20°C (Section 12).

## 12. Storage of the 255X Ph-FI Stock Concentrate for Future Use

If not all of the 255X Ph-FI stock concentrate will be used the same time it is reconstituted, the unused portion may be stored at -20°C for 6 months. During that time, the 255X Ph-FI stock concentrate may be thawed and used twice. After the second thaw, discard any remaining 255X Ph-FI stock concentrate. If

you anticipate using it more than twice, make small aliquots in amber vials or polypropylene tubes and store at -20°C protected from light. Following this protocol, each 500  $\mu\text{L}$  cell sample uses 2  $\mu\text{L}$  of 255X stock concentrate, yielding a 20  $\mu\text{M}$  solution of Ph-FI at 1X. When ready to use, follow Section 13 below.

## 13. Preparation of the 51X Ph-FI Working Solution for Immediate Use

Using the reconstituted 255X Ph-FI stock, prepare the 51X working-strength Ph-FI solution by diluting the stock 1:5 in PBS at pH 7.4. Following the suggested protocols here, each 500  $\mu\text{L}$  sample to be tested requires only 10  $\mu\text{L}$  of 51X Ph-FI solution (or 2  $\mu\text{L}$  of the 255X Ph-FI stock).

1. If you are using the entire vial, simply add 200  $\mu\text{L}$  PBS pH 7.4 to each reconstituted vial (each vial contains 50  $\mu\text{L}$  of the 255X Ph-FI stock; this yields 250  $\mu\text{L}$  of a 51X Ph-FI working solution).
2. If not using the entire vial, dilute the 255X Ph-FI stock 1:5 in PBS, pH 7.4. For example, add 10  $\mu\text{L}$  of the 255X Ph-FI stock to 40  $\mu\text{L}$  PBS (this yields 50  $\mu\text{L}$  of a 51X Ph-FI working solution). Store the unused 255X Ph-FI stock at -20°C (Section 12).
3. Mix by inverting or vortexing the vial at RT.

- **The 51X Ph-FI working solution must be used the same day that it is prepared.**

## 14. General Procedure

Staining cells with ICT's Cholinesterase Detection Kit can be completed within a few hours. However, the kit is used with living cells, which require periodic maintenance and cultivation several days in advance. In addition, once the proper number of cells has been cultivated, time must be allotted for the experimental process. Therefore, as the 51X Ph-FI working solution must be used immediately, it should be prepared at the end of the experimental process. Cells may be manipulated in test tubes, or in plates. The following is a quick overview of the labeling protocol:

1. Expose your cell population to your experimental system (such as induction of apoptosis). Also prepare non-treated negative controls for reference (such as a non-induced cell population).
- **Cell density in the cell culture flasks should not exceed 10<sup>6</sup> cells/mL. Cells cultivated in excess of this concentration may begin to naturally enter apoptosis. Optimal cell concentration will vary depending on the cell line used and your experimental conditions.**
2. Reconstitute each vial of Ph-FI lyophilized reagent (catalog #6156) with 50  $\mu\text{L}$  DMSO and mix thoroughly to solubilize all of the Ph-FI probe contained

in the vial (Section 11). Allow the DMSO to travel around the base of the amber vial until completely dissolved. At room temperature (RT), this reagent should be dissolved within a few minutes. This yields a 255X stock concentrate that may be stored at -20 °C for future use (Section 12), or further diluted to stain the cells (Section 13).

3. Add 200  $\mu$ L PBS to each reconstituted vial and mix (Section 13). This yields 250  $\mu$ L of a 51X Ph-FI working solution that will be used to label the cells. Each 500  $\mu$ L cell sample will use 10  $\mu$ L of the working solution.
4. Label each cell sample using a 1:51 v/v ratio of the working solution. For example, add 10  $\mu$ L of the 51X Ph-FI working solution to 500  $\mu$ L of your cell suspension. This ratio of Ph-FI probe/sample yields a 20  $\mu$ M Ph-FI probe concentration in the labeled sample. You may use a different volume of cells, and your particular cell line may require a different concentration of the Ph-FI probe.

● **When using this assay for the first time, optimize the incubation period for your particular cell line and experimental conditions by setting up several samples to incubate for different lengths of time.**

5. Incubate the labeled cells for 60 minutes at 37°C in a CO<sub>2</sub> incubator protected from light. Resuspend the cells at least once during this incubation period to facilitate even equilibration of the probe into the cells. Active cholinesterase will begin to react with the Ph-FI probe within 15 minutes of addition to the sample. Your particular cell line and experimental conditions may require a longer incubation period (samples typically incubate for 1-4 hours).

● **When using this assay for the first time, set up several samples to incubate for different lengths of time to optimize the incubation period for your particular cell line and experimental conditions.**

6. Dilute the 10X cellular wash buffer 1:10 in DI H<sub>2</sub>O. For example, add 10 mL 10X wash buffer to 90 mL DI H<sub>2</sub>O to make 100 mL (Section 9).
7. After staining the cells and incubating them with the Ph-FI reagent, the cells need to be washed to remove any unbound Ph-FI reagent. The Ph-FI reagent is always fluorescent, so any unbound reagent left with the cells may lead to high backgrounds or false-positives. Wash the cells using ICT's 1X cellular wash buffer, or your cell culture media, or a neutral-pH isotonic PBS-based buffer containing 1% BSA. Cells are washed by removing the labeled media, rinsing the cells in fresh buffer and centrifuging them, or letting them incubate further in fresh buffer (see Sections 15 and 16 below).

● **Several different wash techniques may be used to wash your cells. See Sections 15 and 16 below.**

8. After washing, resuspend the cells in fresh buffer.
9. If not analyzing the cells immediately, cells may be fixed for viewing up to 24 hours later. Add the fixative at a 1:10 ratio into the volume of cell suspension to be fixed. For example, if the cells were resuspended in 500  $\mu$ L, add 50  $\mu$ L fixative.
10. Evaluate the cells. The Ph-FI reagent has an optimal excitation range from 488 - 492 nm, and emission range from 515 - 535 nm (the excitation / emission pairs which best approximate this optimal range should be used).
  - a. Viewing cells through a fluorescence microscope, cholinesterase-positive cells will fluoresce green, while negative cells will appear mostly unstained.
  - b. When analyzing on a fluorescence plate reader (with **black** microtiter plates), use 488 nm excitation coupled with 530 nm emission pairing along with a 515 nm lower wavelength cut-off filter. Cholinesterase activity can be quantitated as the amount of green fluorescence emitted from bound Ph-FI probes. Cell populations with more cholinesterase activity will have a higher RFU intensity than cell populations with less activity.
  - c. Using a flow cytometer, analysis is done using a 15 mW argon ion laser at 488 nm. Fluorescein is measured on the FL1 channel, and a histogram may be generated using the log FL1 (X-axis) versus the number of cells (Y-axis). On this histogram, there will appear two cell populations represented by two peaks. The majority of the negative cells will occur within the first log decade of the FL1 (X) axis (first peak), whereas the cholinesterase-positive cell population will appear as a separate peak or as a shoulder of the first peak showing increased fluorescence intensity.

### 15. Washing Non-Adherent Cell Populations

Non-adherent cells may easily be lost in the wash process, so care is needed when aspirating the media. If growing suspension cells in a plate, the entire plate may be gently spun in the wash process (however this may be harmful to some cell lines so you may want to try several wash techniques). Here is one method of washing the cells:

6. After staining and incubating the cells in Section 14, Step 5, add an extra 1-2 mL of the wash buffer (or other buffer as mentioned above) to each tube.
7. Mix the cell suspension thoroughly.
8. Pellet down the cells by centrifuging them at <300 x g for 5-6 minutes at RT.
9. Carefully remove and discard the supernatant.

10. Resuspend the cell pellets in about 1 mL of wash buffer.
11. Pellet the cells down again by centrifuging them at  $<300 \times g$ , but reduce the centrifugation time down to 3-4 minutes.
12. Repeat steps 9-11.
13. Carefully remove and discard the supernatant.
14. Resuspend the cells in 500  $\mu\text{L}$  of the wash buffer or just PBS or cell media.
15. Fix the cells, if desired (Section 14, Step 9).
16. Evaluate the cells (Section 14, Step 10).

## 16. Washing Adherent Cell Populations

Adherent cells need to be carefully washed to avoid the loss of any cells which round up and come off the plate surface. Loose cells may be harvested from the plate or slide surface and treated as suspension cells, while those remaining adherent to the surface should be washed as adherent cells. The washed non-adherent cells can then be recombined with the adherent portion when the analysis is performed. If growing adherent cells in a plate, the entire plate may be gently spun as part of the wash process (however this may be harmful to some cell lines so you may want to try several wash techniques). Here is one method of washing the cells:

6. After staining and incubating the cells in Section 14, Step 5, pour off or aspirate the Ph-FI-labeled media.
7. Add 1-2 mL of the wash buffer (or other buffer, see Section 14, Step 7).
8. Let the wash buffer incubate with the cells for several minutes at RT protected from light (15 minutes to 1 hour, cells may be placed in the incubator).
9. Harvest the loose cells and wash as suspension cells.
10. Carefully aspirate the supernatant and discard.
11. Wash the cells up to 3 times by repeating Steps 7-9. Some cells may need to be washed more than others depending on the type of instrument used. Cells evaluated by flow cytometry may not need to be washed as much as cells evaluated in a plate reader or microscope, as the sheathing fluid acts as a wash buffer. Cells analyzed in a plate reader often have to be washed the most, as the fluorometer measures total fluorescence within the well and any excess reagent will lead to higher RFUs.
  - a. Cells to be analyzed in a flow cytometer can be trypsinized off the surface of the flask or slide and subsequently resuspended in buffer.
12. Recombine loose cells with adherent cells.
13. Fix the cells if desired (Section 14, Step 9).
14. Evaluate the cells (Section 14, Step 10).

## 17. Sample Microscope Data

**Figure 2**

**Figure 2:** Varying expressions of cellular DNA (A), cholinesterase (B), and caspase activity (C) in apoptotic Jurkat cells. Jurkat cells were treated with 20  $\mu\text{M}$  Ph-FI (B), 10  $\mu\text{M}$  SR-VAD-FMK (C) for 1 hour, and counterstained with 1  $\mu\text{g}/\text{mL}$  of DAPI (A). The cells were viewed using UV- (A), blue- (B) or green- (C) incident light illumination to excite DAPI, Ph-FI or SR-VAD-FMK, respectively. Cell #1 has active cholinesterase (it is green, B), but has little to no caspase activity (it is not red in C), nor does it reveal cellular DNA (it is not blue in A). Cells #2 and #3 reveal cellular DNA and cholinesterase and caspase activity at varied intensities.

### Cholinesterase Detection Kit Ordering Information:

Kit	Reagent	25-test kit	100-test kit
Cholinesterase	Ph-Fluorescein	catalog# 973	catalog# 974

**FLICA™ Ordering Information:** Use the cholinesterase kit with ICT's red fluorescent FLICA™ Apoptosis Detection Kits for dual-labeling studies. The Poly-Caspases kit (#916 and #917) was used in Figure 2:

Caspase	FLICA™ Reagent	25-test kit	100-test kit
Poly-Caspases	SR-VAD-FMK	catalog# 916	catalog# 917
Caspases 3&7	SR-DEVD-FMK	catalog# 931	catalog# 932
Caspase 9	SR-LEHD-FMK	catalog# 960	catalog# 961

## 18. References

1. Park, S.E., N.D. Kim, and Y.H. Yoo. 2004. Acetylcholinesterase plays a pivotal role in apoptosome formation. *Cancer Res.* 64:2652-2655.
2. Zhang, X.J., L. Yang, Q. Zhao, J.P. Caen, H.Y. He, Q.H. Jin, L.H. Guo, M. Alemany, L.Y. Zhang, and Y.F. Shi. 2002. Induction of acetylcholinesterase expression during apoptosis in various cell types. *Cell Death and Differ.* 9: 790-800.
3. Atack, J.R., QS, Yu, T.T. Soncrant, A. Brossi, and S.I. Rapoport. 1989. Comparative inhibitory effects of various physostigmine analogs against acetyl- and butyrylcholinesterases. *J. Pharmacol. Exp. Ther.* 249(1): 194-202.
4. Herholz, K., B. Bauer, K. Wienhard, L. Kracht, R. Mielke, M.O. Lenz, T. Strotmann, and W.D. Heiss. 2000. In-vivo measurements of regional acetylcholine esterase activity in degenerative dementia: comparison with blood flow and glucose metabolism. *J. Neural. Transm.* 107(12): 1457-1468.
5. Giacobini, E. 2004. Cholinesterase inhibitors: new roles and therapeutic alternatives. *Pharmacol. Res.* 50(4): 433-440.
6. Giacobini, E. 2003. Cholinergic function and Alzheimer's disease. *Int. J. Geriatr. Psychiatry* 18(1): S1-5.
7. Ballard, C.G. 2002. Advances in the treatment of Alzheimer's disease: benefits of dual cholinesterase inhibition. *Eur. Neurol.* 47: 64-70.
8. Yu, QS., X. Zhu, H.W. Holloway, N.F. Whittaker, A. Brossi, and N.H. Greig. 2002. Anticholinesterase Activity of Compounds related to geneserine tautomers, N-oxides and 1,2-oxazines. *J. Med. Chem.* 45:3684-3691.
9. Huang, X. B. Lee, G. Johnson, J. Naleway, A. Guzikowski, W. Dai, and Z. Darzynkiewicz. 2005. Novel assay utilizing fluorochrome-tagged physostigmine (Ph-F) to *In Situ* detect active acetylcholinesterase (AChE) induced during apoptosis. *Cell Cycle* 4(1): In Press.